

MEET THE ELECTROLYTES



Electrolyte
Marker

INTRODUCING DIAZYME'S ENZYMATIC ELECTROLYTE ASSAYS

Sodium | Potassium | CO₂

- An alternative to high priced ISE methods
 - Eliminates the need for electrode maintenance
 - Liquid stable reagent and calibrator requires no preparation
 - Accurate and precise
 - Traceable to NIST standard and ISE method
 - Wide range of instrument parameters available for simplifying implementation
 - Low cost per test
- 

ENZYMATIC ELECTROLYTE ASSAYS

Dual Vial
Liquid Stable

	Sodium	Potassium	CO ₂
Method	Enzymatic - A sodium dependent β -galactosidase cleaves ONPG substrate. The product O-nitrophenyl is read at 405 nm and is proportional to the sodium concentration	Enzymatic - A potassium dependent pyruvate kinase catalyzes the conversion of NADH analog to NAD analog which is measured at 380 nm and is proportional to the potassium concentration	The assay is based on a coupled reagent with phosphoenolpyruvate carboxylase (PEPC) and malate dehydrogenase (MDH). The decrease in absorbance at 405 or 415 nm is directly proportional to CO ₂ concentration in the sample
Traceability	NIST standard and ISE method	NIST standard and ISE method	NIST standard and ISE method
Method Correlation to Predicate	$R^2 = 0.9801$ regression $y = 1.0516x - 2.2343$	$R^2 = 0.98$ regression $y = 1.07x - 0.30$	$R^2 = 0.99$ regression $y = 1.0447x - 0.9742$
Precision	The Within-Run CVs 1.2% at 137 mmol/L 1.1% at 160 mmol/L The total CVs 1.56% at 137 mmol/L 1.65% at 160 mmol/L	The Within-Run CVs 1.12% at 4.46 mmol/L 1.20% at 6.86 mmol/L The total CVs 1.77% at 4.46 mmol/L 1.77% at 6.86 mmol/L	The Within-Run CVs 2.3% at 25 mmol/L 2.3% at 40 mmol/L The total CVs 2.8% at 25 mmol/L 3.3% at 40 mmol/L
Reagent On-Board Stability*	4 weeks	2 weeks	4 weeks
Calibrator	Liquid stable calibrator set, no serial dilutions are required	Liquid stable calibrator set, no serial dilutions are required	Liquid stable calibrator set, no serial dilutions are required
Sample Type	Serum	Serum	Serum, Heparinized Plasma
Sample Volume	8 μ L	5 μ L	3 μ L
Linearity	80 - 180 mmol/L	2 - 8 mmol/L	1.12 - 50 mmol/L
Regulatory Status	• 510(k) Cleared •  2797	• 510(k) Cleared •  2797	• 510(k) Cleared • 

*Analyzer Dependent

DIAZYME LABORATORIES, INC.

12889 Gregg Court, Poway, CA 92064 USA
PO Box 85608, San Diego, CA 92186 USA
Tel: +1-858-455-4768 +1-888-DIAZYME
www.diazyme.com sales@diazyme.com

SHANGHAI DIAZYME CO., LTD.

Building 8, 3879 Dongchuan Road
Minhang District, Shanghai, 200245
Tel: +86-21-54843802 Fax: +86-21-51320663
www.lanyuanbio.com service@lanyuanbio.com

DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Germany
Tel: +49-351-886-3300 Fax: +49-351-886-3366
sales@diazyme.de

