

HUMAN KAPPA (κ) AND LAMBDA (λ) FREE LIGHT CHAIN ASSAYS



Cancer
Marker

Diazyme's Human Kappa (κ) and Lambda (λ) Free Light Chain assays provide accurate and reliable testing across the full analytical range and correlates well with current methods. Both assays display excellent precision on the low end while its liquid stable format requires no reagent preparation, saving time and reducing sample handling.

DIAZYME HUMAN KAPPA (κ) AND LAMBDA (λ) FREE LIGHT CHAIN ASSAYS ADVANTAGES

- Assay utilizes polyclonal antibodies for comprehensive detection
- Aids in the diagnosis and monitoring of multiple myeloma in conjunction with other laboratory findings
- Eliminates the need for standalone systems with a wide range of instrument parameters available for simplifying implementation
- Excellent lot to lot performance with a large dynamic range
- No sample pre-dilution is required
- Reduced manual dilutions and retesting for extension modes
- Fast test results (10 minutes) for a rapid turnaround time





REGULATORY STATUS

510(k) Cleared; EU:  ;

Health Canada Registered

HUMAN KAPPA (κ) AND LAMBDA (λ) FREE LIGHT CHAIN ASSAYS

Dual Vial
Liquid Stable



Assay	Kappa (κ)		Lambda (λ)
Method	Latex Enhanced Immunoturbidimetric Assay		Latex Enhanced Immunoturbidimetric Assay
Sample Type & Volume	• Serum, Sample Volume 4 µL		• Serum, Sample Volume 20 µL
Method Correlation	Linear Regression N = 126 y-intercept = -2.536 Slope = 0.958 R ² = 0.977		Linear Regression N = 126 y-intercept = 0.636 Slope = 1.054 R ² = 0.972
Linearity	DZ169A-K (AMI) Standard Mode: 2.9 - 150 mg/L Extended Mode: 2.9 - 3000 mg/L DZ169A-KY1 (AMI)* Standard Mode: 2.3 - 150 mg/L Extended Mode: 2.3 - 3000 mg/L		(AMI) Standard Mode: 3.5 - 200 mg/L Extended Mode: 3.5 - 4000 mg/L
LOB	DZ169A-K	DZ169A-KY1*	1.6 mg/L
LOD	0.9 mg/L	0.9 mg/L	2.8 mg/L
LOQ	1.8 mg/L	1.4 mg/L	3.5 mg/L
	2.9 mg/L	2.3 mg/L	

DIAZYME LABORATORIES, INC.

12889 Gregg Court, Poway, CA 92064 USA
PO Box 85608, San Diego, CA 92186 USA
Tel: +1-858-455-4768 +1-888-DIAZYME
www.diazyme.com sales@diazyme.com

DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Germany
Tel: +49-351-886-3300 Fax: +49-351-886-3366
sales@diazyme.de

SHANGHAI DIAZYME CO., LTD.

Building 8, 3879 Dongchuan Road
Minhang District, Shanghai, 200245
Tel: +86-21-54843802 Fax: +86-21-51320663
www.lanyuanbio.com service@lanyuanbio.com



*Product for Export Only. Not for Distribution in USA.