

HOMOCYSTEINE ENZYMATIC ASSAY



Cardiovascular
Marker



Diazyme's Homocysteine Enzymatic Assay features convenient ready to use reagent, calibrators and controls for the quantitative determination of total L-homocysteine in serum or plasma.

Diazyme's proprietary Enzyme Cycling methodology is an excellent choice for cost conscious laboratories of all sizes due to a wide variety of instrument specific packaging options. The assay required minimal patient sample volume and provides fast, accurate and precise results.

DIAZYME HOMOCYSTEINE ASSAY ADVANTAGES

- Award winning Homocysteine enzymatic assay recognized by the American Association of Clinical Chemistry (AACC)
- Innovative enzyme cycling based technology for accurate and reliable results
- Excellent correlation to HPLC and immunochemical methods
- No "carry over" issues with iron or lipase reagents
- Test renal patients with confidence since there is no interference from cystathionine which affects some other less specific methods
- Available in 2 or 3 reagent format
- Liquid stable format requires no reagent preparation
- Wide range of instrument parameters available for simplifying implementation

REGULATORY STATUS

510(k) Cleared; EU:  ;
Health Canada Registered

HOMOCYSTEINE ENZYMATIC ASSAY

Dual Vial
Liquid Stable

ASSAY SPECIFICATIONS

| | |
|-----------------------------------|---|
| Method | Diazyme Patented Enzyme Cycling |
| Sample Type & Volume | <ul style="list-style-type: none"> Serum Plasma - EDTA - Li-heparin Sample Volume 13 µL |
| Method Comparison | N = 40 y-intercept = 1.05 Slope = 0.94 R ² = 0.99 |
| Linearity | Up to 50 µmol/L |
| Calibration Levels | 5-Point Calibration |
| Prozone/Hook Tolerance | up to 10,000 µg/g |
| Reagent On-Board Stability | Opened: At least 60 days when stored at 2-8°C |

Precision studies were tested with HCY
Enzymatic Assay on OLYMPUS AU400

*Analyzer Dependent

Parameter questions for Homocysteine
Enzymatic Assay should be addressed
to Diazyme technical support. Please call
858.455.4768 or email support@diazyme.com

1. Vilaseca et al. Clin. Chem. 43: 690-692 (1997)

2. Faure-Delanef et al. Am. J. Hum. Genet. 60: 999-1001 (1997)

ASSAY PRECISION

Precision studies were conducted according to the NCCLS
EP-5 protocol. Four HCY serum samples containing 7.0, 12.0,
15.6 and 29.0 µM HCY were tested.

| HCY Concentration | 7 µM | 12 µM | 15.6 µM | 29 µM |
|--|------|-------|---------|-------|
| Within-Run Imprecision CV% N = 20 | 4.5 | 1.87 | 3.04 | 2.4 |
| Total Imprecision CV% N = 30 | 5.87 | 4.88 | 5.51 | 2.57 |

ASSAY INTERFERENCE

An interference study was performed by testing a serum
sample spiked with varied concentrations of endogenous
substances. The following substances normally present in the
serum produced less than 10% deviation when tested at the
stated concentrations:

| | |
|----------------------|------------|
| Bilirubin: | 40 mg/dL |
| Triglycerides: | 1000 mg/dL |
| Hemoglobin: | 500 mg/dL |
| Bilirubin Conjugate: | 40 mg/dL |
| Ascorbic Acid: | 10 mM |
| Cystathionine: | 100 µM** |

**The concentrations tested are about 5-10 times higher than
the normal range of serum levels.

REFERENCE RANGE

In most of the U.S. clinical laboratories, 15 µmol/L is used
as the cut-off value for normal level of Hcy for adults.¹⁻² In
Europe, 12 µmol/L is used as the cut-off value. However, each
laboratory is recommended to establish a range of normal
values for the population in their region.

DIAZYME LABORATORIES, INC.

12889 Gregg Court, Poway, CA 92064 USA
PO Box 85608, San Diego, CA 92186 USA
Tel: +1-858-455-4768 +1-888-DIAZYME
www.diazyme.com sales@diazyme.com

DIAZYME EUROPE GMBH

Zum Windkanal 21, 01109 Dresden, Germany
Tel: +49-351-886-3300 Fax: +49-351-886-3366
sales@diazyme.de

SHANGHAI DIAZYME CO., LTD.

Building 8, 3879 Dongchuan Road
Minhang District, Shanghai, 200245
Tel: +86-21-54843802 Fax: +86-21-51320663
www.lanyuanbio.com service@lanyuanbio.com

